



BETAJECT® Comfort

How to get your BETAJECT® Comfort

To order the BETAJECT® Comfort, please contact your local Schering organisation or Betaferon® nurse.

For further information on Betaferon®, please refer to our MS Information Gateway website at www.ms-gateway.com.




EU2006.0817



The new generation of autoinjectors

Suitable for the new Betaferon® injection system. Autoinjectors have been proven to reduce injection site reactions. They are particularly suitable for people with a fear of needles.



The new BETAJECT® Comfort

Customised

New three-step injection technology based on the natural injection process:

- 1 Needle insertion into the skin
- 2 Slow release of the medication
- 3 After completed injection, gentle needle extraction back into the autoinjector

Flexible

- Flexible depth adjustment possible
Please ask your health care professional for advice
- Suitable for dose deliveries from 0.25 to 1.0 ml
- Can be used with the thinnest needle size available (30 G)

Convenient


- A bell as well as a new visual marker indicating end of the completed injection process
- Prolonged injection time for gentle administration
- Comes in an attractive case

Ready, Steady, Go

Step 1: Load the device



Step 2: Insert the ready-to-use syringe



Step 3: Inject



The benchmark autoinjector with the three-step technology.